



THE ESSENTIAL GUIDE TO  
**Selecting Your PT Electrode Tips**  
**Manufacturing Partner**



## Table of Contents

Selecting Your PT Electrode Tips Manufacturing Partner	3
Strategic Investment & Supply Chain Optimization	4
Importance and Versatility of Electrophysiology Catheters	5
Manufacturing Challenges of Electrode Tips for Electrophysiology Catheters	6
Material Expertise in Platinum-Iridium	6
Precision Machining Capabilities	7
Stringent Quality Control	8
Risk Management and Metrology	8
Challenges in Achieving Uniform Surface Finish	9
Pulse Technologies Joins Integer	10
Pulse Technologies at a Glance	11

# Selecting Your PT Electrode Tips Manufacturing Partner

The global market for electrophysiology catheters valued at \$8 billion in 2024 is poised for rapid expansion, projected to grow at high single digit growth through 2030. This expansion is propelled by several key factors: the rising incidence of cardiovascular disorders such as arrhythmias and atrial fibrillations, an increasing shift towards minimally invasive procedures, and significant technological advancements. These technological strides include the development of 3D mapping systems and pulsed field ablation, both of which are revolutionizing the treatment of heart arrhythmias.

1

### Increasing Prevalence of Cardiovascular Disorders:

Cardiovascular diseases remain a leading cause of death globally, with arrhythmias and atrial fibrillations being particularly prevalent. The growing number of patients diagnosed with these conditions is a primary driver of demand for advanced diagnostic and therapeutic tools in electrophysiology.

2

### Preference for Minimally Invasive Procedures:

There is a strong trend toward minimally invasive procedures, which reduce recovery times and improve patient outcomes. These procedures, including those using electrophysiology catheters, are becoming the standard of care for treating heart arrhythmias.

3

### Technological Advances in Electrophysiology:

Advances such as 3D mapping systems provide clinicians with detailed and accurate images of the heart's electrical activity, enabling precise targeting during ablation procedures. Pulsed field ablation (PFA) represents a breakthrough in how cardiac arrhythmias are treated, offering a non-thermal method that minimizes tissue damage and enhances patient safety.



## Strategic Investment & Supply Chain Optimization



The electrophysiology market's substantial growth is fueled not only by technological advancements but also by robust funding for research and proactive government initiatives aimed at enhancing patient care. These efforts are pivotal in fostering innovations and improving treatment methodologies in electrophysiology, thereby contributing significantly to the market's expansion.

For medical device companies navigating this rapidly evolving landscape, strategic partnerships with adept contract manufacturers are indispensable. Such collaborations are essential to scaling operations efficiently, meeting the surging demand, and elevating customer service levels. A responsive and flexible supply chain is critical, particularly one that focuses on components like electrode tips, which are central to the functionality of electrophysiology catheters.

## Importance and Versatility of Electrophysiology Catheters


Electrophysiology catheters are indispensable tools in modern cardiac care, serving both diagnostic and therapeutic functions with high precision and reliability. These catheters are intricately designed to navigate the complex geometry of cardiac anatomy, enabling clinicians to gather essential cardiac data that is critical for accurate diagnosis and targeted treatment.

- **Diagnostic Excellence:** Primarily, electrophysiology catheters are utilized for detailed cardiac mapping. Their ability to adeptly maneuver within the heart allows for the precise capture of vital electrophysiological signals. This capability is crucial for assessing the electrical activity of the heart and pinpointing areas of abnormal conduction. By providing a comprehensive view of cardiac electrophysiology, these catheters facilitate the development of effective treatment plans and enable more precise therapeutic interventions.
- **Therapeutic Applications:** Beyond diagnostics, electrophysiology catheters play a pivotal role in the treatment of cardiac arrhythmias. They are particularly effective in procedures such as catheter ablation, where targeted areas of the heart tissue are modified to correct abnormal electrical signals. This direct intervention helps restore normal rhythm and function to the heart, significantly improving patient outcomes.
- **Versatility in Cardiac Arrhythmia Management:** Electrophysiology catheters also excel in the comprehensive collection and analysis of intracardiac electrograms. This functionality is essential for diagnosing complex arrhythmias and for guiding therapeutic interventions during ablation procedures. Additionally, these catheters can be used for cardiac stimulation and temporary pacing during diagnostic or therapeutic procedures, further underscoring their versatility and critical role in managing heart rhythm disorders.



## Manufacturing Challenges of Electrode Tips for Electrophysiology Catheters

Electrode tips are vital in diagnosing and treating heart conditions — essential for sensing, recording, and stimulating the heart's electrical activity. Precision in the design and manufacturing of these tips significantly impacts the quality of signals they capture and transmit, which in turn affects the accuracy of diagnoses and the effectiveness of treatments. The right contract manufacturer must not only ensure the electrode tips' efficiency and reliability but also their consistent performance, which is crucial for the overall effectiveness of the catheters. As such, choosing a manufacturing partner who can provide high-quality, dependable electrode tips is a key determinant of a catheter's clinical success.



## Material Expertise in Platinum-Iridium

Electrode tips are commonly manufactured from platinum-iridium, a high-performance alloy chosen for its biocompatibility, corrosion resistance, and excellent electrical conductivity. However, working with platinum-iridium presents specific challenges due to its high melting point and considerable hardness. These properties necessitate a deep understanding of the material's behavior under various machining processes to maintain the integrity of the alloy and achieve the desired properties in the final product.



## Precision Machining Capabilities

Manufacturing electrode tips includes the development of intricate geometries and micro-scale features that demand advanced machining capabilities. The hardness of platinum-iridium requires specialized tools and techniques, such as laser machining, Electrical Discharge Machining (EDM), and Computer Numerical Control (CNC) machining centers. These technologies help ensure the precise shape and dimensions are consistently achieved.

Moreover, it is essential that the contract manufacturer has a robust understanding of Design for Manufacturability (DFM) and Design for Assembly (DFA) principles. These principles guide the design and production processes to simplify and optimize the manufacturing without compromising the design's complexity and functionality.



## Stringent Quality Control

Ensuring the dimensional accuracy and overall quality of electrode tips is essential for their performance in medical applications. Given their complex geometries, inspecting and verifying the dimensions and surface finish of each component presents significant challenges. Manufacturers must be equipped with advanced inspection techniques capable of identifying minute defects such as burrs, scratches, pits, and irregularities in irrigation holes.

Key quality control measures include adherence to FDA's 21 CFR 820 and ISO 13485 regulations, a well-documented Quality Policy, and compliance with Good Manufacturing Practices (GMPs). Furthermore, manufacturers should have robust process validation capabilities to ensure manufacturing repeatability and reproducibility, including performing Installation Qualifications (IQ), Performance Qualifications (PQ), and Operational Qualifications (OQ).

## Risk Management and Metrology

Proactive risk management is critical in the manufacturing process. Utilizing tools like Process Failure Mode and Effects Analysis (PFMEA) helps in anticipating potential failures and implementing strategies to mitigate them. Additionally, having an in-house metrology lab equipped with the latest calibration and measurement tools allows for continuous monitoring and adjustment of production processes to maintain high-quality standards.

## Challenges in Achieving Uniform Surface Finish

The surface finish of electrode tips is crucial not only for their performance but also for their biocompatibility and electrical efficiency. Achieving a uniform metal plating surface that adheres tightly to the base material without peeling or unevenness is critical for the reliability and safety of medical devices. The complex and specialized geometry of these components often complicates the plating process, requiring high levels of expertise and precise control over the process conditions.





## Pulse Technologies Joins Integer

Pulse Technologies has recently joined forces with Integer, one of the world's largest medical device contract development and manufacturing organizations (CDMO), serving diverse markets including cardiac, neuromodulation, vascular, and portable medical devices. This strategic alignment significantly enhances reach and capabilities, bringing several key benefits and capabilities to clients, including:

- **Single-Source Manufacturing Solution:** Streamlining the production process by offering a comprehensive range of manufacturing services under one roof.
- **Large Capacity:** Ability to scale production quickly to meet high demand without sacrificing quality.
- **Dual Sourcing Capability:** Enhances supply chain security by providing multiple sourcing options, reducing the risk of disruption.
- **Integrated Manufacturing:** Coordination from design to production to assembly, ensuring seamless project management and consistency.
- **Rapid Response Times:** Agile and responsive manufacturing processes that can adapt to sudden changes in market demand or project scope.
- **Customization Capability:** Advanced technical skills and equipment to customize products according to specific customer requirements.
- **Low-Volume High Mix Expertise:** Specialized in managing and producing a diverse range of products in smaller quantities to meet unique client needs.
- **Supply Chain Resilience:** Robust supply chain strategies are in place to ensure reliable delivery schedules and minimize supply-related uncertainties.



## Pulse Technologies at a Glance

With over 30 years of experience in medical device manufacturing, Pulse Technologies offers unmatched expertise and a proven track record of innovation and quality:

- **Industry-Leading Metrology and Testing Capabilities:** Our in-house metrology lab is equipped with advanced testing and measurement technologies, ensuring products meet the highest standards.
- **70,000-Square-Foot Facility:** Ample space equipped with state-of-the-art machinery to handle projects of any scale.
- **100+ CNC Machine Tools:** A wide array of sophisticated CNC machinery that allows for precision manufacturing of complex components.
- **Data Management:** Expertise in formatting data for Design History Files (DHF) and Device History Records (DHR), streamlining documentation and compliance.
- **For more information about our comprehensive manufacturing capabilities and how we can support your medical device production needs, please contact us today.**



**Integer | Corporate Headquarters**

5830 Granite Parkway, Suite 1150  
Plano, TX 75024  
contact-integer@integer.net  
[Integer.net](http://Integer.net)

**Manufactured by Pulse Technologies, Inc. – An Integer Company**

P: 267-733-0200 | F: 267-733-0201  
[sales@pulsetechnologies.com](mailto:sales@pulsetechnologies.com)  
[www.pulsetechnologies.com](http://www.pulsetechnologies.com)